
*Bob
AND
Brad*[®]

User Manual

QL/DMS.C2-D
Mini Massage Gun

Maintenance

The battery in this device is not designed to be removed or exchanged by the user. If any problem with the power supply occurs, please contact us directly.

Do not try to open the device, repair it or fix it by taking the device apart yourself. Tampering with the device's mechanics in any way will result in an automatic expiration of the product 's warranty.

Warnings and Safety Concerns

Before connecting the device to the power supply, please ensure that the supply voltage stated on the rating plate is compatible with the unit. Please charge the new batteries before the first use.

Do not submerge the device in water or any other liquid.

Do not use the device for more than 30 minutes at a time. Wait 5 minutes prior to using the device again, following a 30-minute session with the device.

Do not apply to your head or near your genitals.

Do not use the device to supplement or replace any medical treatment. Please consult your physician before using the device if you are experiencing unexplained pain, if you're receiving medical treatment, or using medical equipment.

FCC Statement

1. This device complies with Part 15 of the FCC Rules. Operation is subject to the following

two conditions:

- (1). This device may not cause harmful interference.
- (2). This device must accept any interference received, including interference that may cause undesired operation.

FCC Statement NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These

limits are designed to provide reasonable protection against harmful interference in residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause any harmful interference to radio or television receptions, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help with the supplied fastener. **DO NOT** hang the unit by the power cord.

DISCLAIMER


Please contact your physician before using the Bob and Brad **QL/DMS.C2-D** if you have any physical health concerns or have any of the following conditions:

pregnancy, diabetes, epilepsy, migraines, herniated disks, spondylolisthesis, recent joint replacements or IUDs. Please also contact your physician before use if you have undergone recent surgery or have any metal pins or plates in your body.

Contact us

www.bobandbrad.com
support@bobandbrad.com
Tel:+1 (612)567-3035
M-F 9:00AM-5:30PM

Main structure and the fittings

Operation steps

1. Connect the device with the massage head when the device is turned off, insert the massage head completely with pressure and pull it out when removing the vibrating head (note: The Fork head should be positioned for installation. When the buckle is inserted into the slot, it means the installation is in place. If the massager head is not installed in place, there will be a risk of falling while in use.)

2. Turning the device ON/OFF: Press and hold the " ⏻ " power button for 2 seconds.
3. To change the percussion massage frequency - when the device is on, press the button " ⏻ " to turn on the first setting. If you would like to advance to a higher setting, press the power button " ⏻ " again. The active setting will be indicated by the back light. To pause the massage gun press the power button for 2 seconds. Press for 2 seconds again to resume.

Note:

- A. After 10 minutes of operation, the device automatically turns off.
- B. LED light indicates the corresponding strike frequency level, the lowest at level 1 and the highest at level 5.
- C. When the device is standby mode (0 level), it will automatically shut down without any operation for 20 seconds.

Battery display

| LED1 | LED2 | LED3 | LED4 |
|------------|------------|------------|-------------|
| ● | ● | ● | ● |
| 25% | 50% | 75% | 100% |

***LED1 flashing means the battery is running down to less than 10%**

When the first indicator light is on, it's indicating 25% of the battery is remaining.

When the second indicator light is on, it's indicating 50% of the battery is remaining.

When the third indicator light is on, it's indicating 70% of the battery is remaining.

When the fourth indicator light is on, it's indicating 100% of the battery is remaining.

Note:

1. When the first indicator light flashes the battery is running down to less than 10% and a recharge is required.
2. It is not recommended that the battery power be completely depleted to 25%.

Charging

1. Before the first use, please fully charge the device and to charge fully it takes about 3 hours.
2. Connect the output end of the connection line to the DC charging port at the bottom of the device, and then connect the input end of the connection line to the mobile power source or mobile phone charger (5V/2A, if the charging current or voltage is less than 5V/2 A, charging time will be longer. Any regular phone adapters will be able to charge the massage gun without a problem.)
3. The battery icon on the handle flickers regularly. The digital display on the handle shows the battery power and indicates when the battery is being charged and the level of power remaining.
4. When the battery icon on the handle stops flickering and remains on, the digital display shows that 100% of the battery is full. The battery can be recharged at any level of battery power. When charging, please ensure that the device is shutdown.

Technical details:

QL/DMS.C2-D Mini Massage Gun

Speed: 2000 – 3200 r/min

Rated Power: 25W

Dimensions: 18.6x14.7x6cm / 7.3" x5.7" x2.3 "

Battery Information:

Rated Voltage: 7.4V 

Battery capacity: 2500mAh

Charging time: 210 minutes

Guarantee

This product comes with a one-year warranty from the date of purchase.

Warranty will be voided due to any of the following:

Damaged due to improper usage or storage.

Damaged by attempting to repair or disassemble the device. Period of warranty exceeded.


Utilize the Bob and Brad massage gun before, during, or after your workout for maximum benefits. This massage gun can help accelerate muscle tissue responsiveness by increasing hydration, range of motion and function to tight muscles.

| Indication | Action | Speed | Pressure | Duration | Muscle Position | Movement |
|-----------------------------------|--|--------------|--|-----------------|---------------------------------|--|
| Warm Up / Activation | Activate and prime muscles prior to exercise. Reduce risk of injury. | Level 3 | Light to Moderate | 30 sec or less | Relaxed. Move along the muscle. | Quick, high speed over area. |
| Mid Workout / Myofascial Mobility | Help muscle recovery and reactivation. | Level 4 | Heavy (especially for large muscle groups) | 8-12 sec | Target specific area. | Localised |
| Recovery | Relax muscles and aid recovery. Down-regulate or turn off muscles. | Level 1 or 2 | Moderate to Heavy | 90 sec or more | Move along the muscle. | Slow movement along the muscle |
| Pain Modulation | Reduction in pain to increase movement. | Level 1 | Light | 10 to 30 sec | Relaxed. Localised spot. | Over painful area. |
| Body Awareness | Improve movement or mechanics of a joint. Activation of a muscle. Draw attention to a muscle | All Levels | Light to Moderate | 10 to 20 sec | Relaxed or with movement | Over specific target area or movement over area. |


Full Body Massage

BACK

The person receiving the back massage should remain relaxed and lie on a comfortable massage table or soft surface on the floor.


- ▷ **Fork Head Attachment:** Put both ends of this attachment along both sides of the spine, vibrate and massage along the whole spinal area. Level 1 vibration speed is recommended. Please avoid placing on the spine directly.

- ▷ **Flat Head Attachment:** Massage the left and right waist gently and repeatedly. Level 2-3 vibration speed is recommended. The effect of this massage gun on the treatment of kidney stone disease has not been verified by clinical trials. Please consult your doctor before use.

- ▷ **Flat or Air Cushion Head Attachment:** For massaging on the scapula and back. Level 1 vibration speed is recommended.

- ▷ **The Ball or Flat Head Attachment:** Massage the muscle tissue on both sides of the lower waist. Level 2-3 vibration speed is recommended.

Full Body Massage


NECK


- ▷ The ball, air cushion and flat head attachments are recommended for neck and shoulder muscles. Please avoid using on the skull. We recommend to use the air cushion attachment for thin or less muscled areas to avoid pain or discomfort.

ARM & PALM

- ▷ The ball, air cushion and flat head attachments are recommended for massaging the arm muscles. For thinner clients, please use the air cushion head on Level 1 vibration speed to avoid pain or discomfort. For overweight or clients with well-developed muscles, using the round head attachment on Level 2-3 speed improves the treatment effect.


- ▷ The bullet head attachment is recommended for treatments of the palm area, especially on small muscle tissues and joint parts. Level 1 - 2 is recommended.

Full Body Massage

LEG & FEET

- ▷ The ball and air cushion attachments are recommended for hip muscles on Levels 2-3 speed to improve the treatment effect. Please avoid using on the spine or tailbone.

- ▷ The ball and air cushion head attachments are recommended for thigh muscles on vibration speed Levels 2-3.

- ▷ Above mentioned head attachments are also applied for a deeper massage on the back thigh and calf muscles.


- ▷ The bullet head attachment is recommended for small muscle tissues on the arch and joint parts. Vibration speed Levels 1 - 2 is therefore recommended.

Full Body Massage

ABDOMEN


During massaging the abdomen, please stay in a standing posture. Level 3 vibration speed is recommended for the abdomen area. The purpose is to use fast shockwaves to vibrate adipose tissue under the skin. Do not use the vibration device on the abdomen continuously for more than 5 minutes to prevent an adverse reaction of the digestive system. Usage on the abdomen is to achieve passive exercise, weight loss, and body shaping treatment. A vibration massage on the abdomen is recommended after sport.


- ▷ The air cushion head attachment is recommended on the abdomen area and adipose tissue on stomach.

Please do not massage your abdomen immediately after dinner.

MORE INFORMATION ABOUT A VIBRATION MASSAGE


- ▷ No matter which head attachments used, shockwaves of this device are able to vibrate through skin and surface tissues. During the treatment, adding extra pressure on the massage device will increase the intensity and depth of shockwaves therefore improving the treatment without influencing the vibration frequency.

- ▷ To increase the intensity and depth of shockwaves, it is recommended to hold the muscle parts around the treated area steadily with hands. By preventing the diffusion of shockwaves through skin and surface tissues, the massage gun can achieve a deeper effect on muscle parts.

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